

Quickclean®

Expert of Pain Relief of OA

SINGCLEAN

How do I get more information
about Quickclean®?

Please email: marketing@singclean.net

Model No.

| Product | Specifications | Model No. --Quickclean® |
|---|------------------|-------------------------|
| Medical Sodium Hyaluronate Gel for Bone Joint | 10 mg/mL: 1.0 mL | ORT-QC-10-1.0 |
| | 10 mg/mL: 2.0 mL | ORT-QC-10-2.0 |
| | 10 mg/mL: 2.5 mL | ORT-QC-10-2.5 |
| | 12 mg/mL: 1.0 mL | ORT-QC-12-1.0 |
| | 12 mg/mL: 2.0 mL | ORT-QC-12-2.0 |
| | 12 mg/mL: 2.5 mL | ORT-QC-12-2.5 |
| | 20 mg/mL: 2.0 mL | ORT-QC-20-2.0 |
| | 20 mg/mL: 3.0 mL | ORT-QC-20-3.0 |

CE
2460



 Quickclean®

Medical Sodium Hyaluronate Gel

Dedicated to Long-term Pain Relief of
OA for 6 Months

Hangzhou Singclean Medical Products Co., Ltd.

Add: No.125(E),10th Street, Hangzhou Qiantang Area, Zhejiang, China

P.c.: 310018

Tel: +86-571-63431868

WhatsApp: +8613362195396

Fax: +86-571-63431886

E-mail: marketing@singclean.net

www.singclean.net


Quickdean®



What Is Osteoarthritis (OA) of the Knee?

OA is a common condition that causes joints to become stiff and painful. It is most common in the elderly.

With OA of the knee, less synovial fluid which contains hyaluronic acid is produced naturally in the knee joint while synovial fluid acts as the knee's shock absorber and protects cartilage coverings from wear and tear. Gradually, as joint fluid continues to break down and cartilage deteriorates, bones may start to rub against each other, leading to pain, stiffness, swelling and limited joint activity.



How Do I Know If I Get OA?

The symptoms of osteoarthritis include

- Joint Pain
- Joint Swelling
- Joint Stiffness
- Reduced Joint Range of Motion
- Joint Warmth
- Audible Sounds
- Deformities

If you are experiencing any or all of the above symptoms, consult with your doctor immediately. The earlier OA is diagnosed, the sooner patients can receive appropriate treatment.



Who Is at Risk for OA?

- People over age 45
- Obesity, overweight will put an extra pressure on knees
- Gender, OA is more prevalent in women, after age 55
- Young people with sports related injuries
- Athletes needing viscoelastic supplement for joint lubrication and mechanical support
- People with osteoarthritis have already been treated with pain relievers and other treatments did not work well

Why Is Hyaluronic Acid Injection Effective for OA?

HA is a natural component of synovial fluid. In OA, the synovial fluid may become less viscous and supply less lubrication to the joints, leading to pain and discomfort.

Hyaluronic Acid Injection can supplement this lost synovial fluid, provide lubrication and cushioning, reduce inflammation, and further protect cartilage and other joint tissues to help alleviate knee pain and restore mobility.

What Is Quickclean®?

A POWERFUL VISCOSUPPLEMENTATION TO HELP YOU RESTORE ENERGY AND YOUTH



Quickclean®, is a kind of hyaluronate from streptococcal fermentation metabolites and is injected directly into your knee for OA treatment by your healthcare provider.

By replacing and supplementing the pathological synovial fluid in the osteoarthritic joints, it reduces pain and improves joint function.

Quickclean® is intended for individuals who find simple pain relievers or other conservative therapies, such as exercise and physical therapy, insufficient for their relief.

Quickclean® Series Details


| Features | Quickclean® | Quickclean® Plus |
|------------------------------------|---|--|
| Picture | | |
| Concentration | 10mg/ml, 12mg/ml | 20mg/ml |
| Viscosity Average Molecular Weight | ≥1.0 Mda | ≥1.6Mda |
| Volume | 1.0ml, 2.0ml, 2.5ml | 2.0ml, 3.0ml |
| Type of HA configuration | Linear | Linear |
| Monthly administration | 3 - 4 injection, once a week for 3 to 4 consecutive weeks | 1 Injection |
| Duration | Up to 6 months | 6 months+ |
| CE | Yes | Yes |
| Who is suitable | For People who like a series of treatments for financial or other reasons | For People who require faster recovery |

Quickclean® Clinical Data

“Quickclean® has been proven to provide long-lasting relief of OA knee pain.”

A total of 300 cases of arthroplogosis subjects were enrolled, the clinical study was performed in First Affiliated Hospital of Zhejiang Chinese Medical University and Tongde Hospital of Zhejiang Province, two Grade A hospitals in China.

Average Effectiveness Rate of Quickclean®



• Performance

In clinical studies, Quickclean® has shown to relieve OA knee pain for 6 months.

The **average effectiveness rate** was 99.33% of 20mg/ml product and 99.00% of 10mg/ml product during 6 months.

• Safety

No adverse events occurred.
No serious adverse events occurred.

Quickclean® Advantages



Better Patient Compliance

Excellent biocompatibility due to bio-fermentation, no immunogen, biodegradable.



Effective

- High Concentration of HA: Creating a high viscosity environment to replicate healthy joint function.
- High Molecular Weight
- Long lasting pain relief
- Delay total knee replacement (TKR) surgery



Convenient

Administered in doctor's office

Quickclean® Popular in the World



Quickclean® International Registration

Quickclean® Medical Sodium Hyaluronate Gel has obtained registration in over **40 countries**, including **UK, European Union, Malaysia, Paraguay, Venezuela, Egypt, Türkiye** and many more.

With its expanded reach, Quickclean® is bringing relief to an even larger population of individuals suffering from OA.

Quickclean® Wonderful Feedback

Positive feedback from patients, which includes reports of improved mobility and pain relief, perfectly underscores Quickclean® 's efficacy.